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

Flow cytometric analysis of non-transfected or His fusion gene transfected CHO cells using **THE™ His Tag Antibody, mAb, Mouse (GenScript, A00186)** (black and green, respectively).

The signal was developed with FITC conjugated Goat Anti-Mouse IgG.

Western blot analysis of cross-reactivity of **MonoRab™ Anti-His Tag (C-term) Antibody (25B6E11), mAb, Rabbit (A01857-40)**

Primary Antibody:

A: MonoRab™ Anti-His Tag Antibody (A01857-40) 1 µg/mL


B: Rabbit IgG (Negative control) 1 µg/mL

Secondary Antibody:


Goat anti-Rabbit IgG (H&L) [IRDye800] (Licor,926-32211)

0.125 µg/mL

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Western blot analysis of C-terminal Strep II tagged fusion protein using **THE™ NWSHPQFEK Tag Antibody, mAb, Mouse (GenScript, A01732, 0.2 µg/mL)**. The signal was developed with IRDye™ 800 Conjugated affinity Purified Goat Anti-Mouse IgG.

Sensitivity and specificity of **MonoRab™ DYKDDDDK Tag**

**Antibody, mAb, Rabbit (A01868-40)** to DYKDDDDK-tag fusion protein by Western blot Loading:

A: N-terminal DYKDDDDK-tagged fusion protein (25 ng, 10 ng, 5 ng)

B: C-terminal DYKDDDDK-tagged fusion protein (25 ng, 10 ng, 5 ng)

C: Internal positions DYKDDDDK-tagged fusion protein (50 ng, 25 ng, 10 ng, 5 ng)


Primary Antibody:

MonoRab™ DYKDDDDK Tag Antibody, mAb, Rabbit (A01868-40) 0.2 µg/mL

Secondary Antibody:


Goat anti-Rabbit IgG (H&L) [IRDye<sup>800</sup>] (Licor,926-32211) 0.125 µg/mL

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
Consistency analysis of Batch 1#, 2#, 3# and 4# of **THE™ His Tag Antibody, mAb, Mouse** (GenScript, **A00186**, 1 µg/mL) by Western blot, showing that signal remains consistent from Lot to Lot. The assay was performed with His-tagged fusion protein. The signal was developed with IRDye™ 800 Conjugated Goat Anti-Mouse IgG.

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Immunocytochemistry/Immunofluorescence analysis of non-transfected CHO cells (a) or His-tagged protein transfected CHO cells (b)

Using **THE™ His Tag Antibody [iFluor 488], mAb, Mouse** (**A01800-100**, 4 µg/mL).

Flow cytometry analysis of non-transfected CHO cells (a) and V5 fusion gene transfected CHO cells (b) using **THE™ V5 Antibody, mAb, Mouse** (**GenScript, A01724**). The signal was developed with FITC conjugated Goat Anti-Mouse IgG.

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产品类型	产品名称	货号
Resin	Ni-IDA Resin	L00223I
	Ni Resin FF	L00465
	High Affinity Ni-Charged Resin FF	L00666
	MonoRab™ Anti-DYKDDDDK Affinity Resin	L00766
	Anti-DYKDDDDK G1 Affinity Resin	L00432
	Anti-DYKDDDDK IP Resin	L00425
	Glutathione Resin	L00206
Magnetic Beads	Ni-charged MagBeads	L00295
	AmMag™ Ni Magnetic Beads	L00776
	MonoRab™ Anti-DYKDDDDK Magnetic Beads	L00835
Kit	GST Fusion Protein Purification Kit	L00207

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产品名称	货号
His Tag ELISA Detection Kit	L00436
His Tag Antibody Plate (black, 96 wells, non-removable)	L00440B
His Tag Antibody Plate (clear, 8 wells × 12 strips)	L00440C
His Tag Antibody Plate (white, 96 wells, non-removable)	L00440W
DYKDDDDK Tag Antibody Plate (black, 96 wells, non-removable)	L00455B
DYKDDDDK Tag Antibody Plate (clear, 8 wells × 12 strips)	L00455C
DYKDDDDK Tag Antibody Plate (white, 96 wells, non-removable)	L00455W
GST tag ELISA Detection Kit	L00411

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产品名称	货号
Enterokinase, His, Bovine	Z03004
TEV Protease, His	Z03030
HRV-3C Protease, His Tag	Z03092

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