

Taq DNA Polymerase without Mg2+

DATASHEET

Version: 201807

Cat. No.: E00008 Size: 1,000 U Description:

Taq DNA Polymerase is a thermostable DNA Polymerase isolated from an *E. coli* strain that carries the *Taq* DNA polymerase gene. *Taq* DNA Polymerase is the most common polymerase used for PCR* reactions.

Key Feature:

Terminal transferase activity. Taq DNA Polymerase has terminal transferase activity which results in the addition of a single nucleotide (adenosine) at 3'end of the extension product.

High-purity. No contamination activity has been detected in standard test reactions.

Unit Definition:

One unit is the amount of enzyme that can incorporate 10 nmol of dNTP into acid-insoluble material in 30 minutes at 74°C.

10 X reaction Buffer (without Mg2+):

500 mM KCl, 100 mM Tris HCl (pH 9.0 at 25°C), 1% Triton X-100 Buffer. This buffer is optimized for use with 200 µM dNTPs.

Note: If the reaction is performed without this buffer, then add 0.1% Triton X-100 (final concentration) to ensure high activity.

Storage Buffer and Concentration:

The enzyme is delivered in 5 units/µl in 20 mM Tris HCl (pH 8.0), 0.1 mM EDTA, 1 mM DTT, 0.1% Triton X-100 and 50% glycerol.

Storage:

This product can be stored at -20°C for future using.

Formulation:

GenScript *Taq* DNA Polymerase has been formulated using GenScript's proprietary technology. The enzyme can be shipped at room temperature or even 37°C for seven days without any loss of activity.

Applications:

Taq DNA Polymerase can be used in most applications including the following:


PCR*

3' A-tailing of blunt ends

Primer extension

DNA sequencing

PCR performance, activity, nuclease.



Lane	Taq	United Used
1		0.1
2	Leading Brand A	0.25
3		0.5
4	Leading Brand B	0.1
5		0.25
6		0.5
7		0.1
8	Genscript	0.25
9		0.5

^{*} The PCR process is covered by US. Patents Nos. 4683195 and 4683202, issued to Cetus and owned by Hoffman-La Roche Inc. GenScript does not encourage or support the unauthorized use of the PCR process. Use of this product is recommended for persons who either have a license to perform PCR or are not required to obtain a license. Sale of this product is restricted to regions or countries where native Taq DNA polymerase patents have been invalidated.