

Version: 03
 Update: 04/28/2022

DATASHEET

DNase I

Cat. No.: E00053-100; E00053-500; E00053-1000

Size: 100 U/ 500 U/ 1000 U/ 5000 U/ 10000 U

Product Introduction

Deoxyribonuclease I (DNase I) is DNA-specific endonuclease that cleaves both single-stranded DNA, double-stranded DNA and DNA-RNA hybrids, yielding 5'-phosphate-terminated polynucleotides with a free hydroxyl group on position 3'-producing tetranucleotides. The activity of DNase I is strictly dependent on Ca^{2+} and can be activated by divalent metal ions such as Mg^{2+} or Mn^{2+} . In the presence of Mg^{2+} , DNase I nonspecifically recognizes and cleaves a double-stranded DNA at any site on either strand, and in the presence of Mn^{2+} , it recognizes and cleaves almost the same sites on both strands of the DNA to produce DNA fragments with blunt ends or sticky ends with 1~2 nucleotide overhangs.

GenScript is offering DNase I produced by expression in a *P. pastoris* strain carrying a plasmid encoding the bovine DNase I.

Source: Recombinant DNase I expressed by yeast.

Species: Bovine

Molecular Weight: 35-40 kDa, on SDS-PAGE under reducing conditions.

Unit Definition: One unit is the amount of the enzyme that increases the absorbance at 260 nm by 0.001 per minute at 25 °C, pH 5.0, with calf thymus DNA as the substrate.

Optimal active temperature: 37 °C

Formulation: Supplied as a solution of 20 mM sodium acetate, 5 mM $CaCl_2$, 0.1 mM PMSF, 50% (v/v) glycerol, pH 6.5 at 25 °C.

Storage & Stability: This product remains stable up to 12 months at -20 °C. Avoid repeated

freeze-thaw cycles.

Inactivation: Add EDTA with a final concentration of 2.5 mM and heat the solution at 65 °C for 10 min can inactivate DNase I.

Application:

- DNA template digestion following *in vitro* transcription
- Genomic DNA digestion prior to RT-PCR
- Preparation of DNA-free RNA samples
- Nick-translation
- Studies of DNA-protein interactions by DNase I, RNase-free footprinting
- Prevent cell clumping without affecting cell viability

Quality Control Specifications

Assay	Specifications
Appearance	Clear, colorless liquid
Purity	≥ 95% as analyzed by SDS-PAGE
Enzyme Activity	≥ 2 U/μl
Endotoxin Level	≤ 0.1 EU/μg of protein by gel clotting method

Residual RNase	Non-detectable
----------------	----------------

Reagents Supplied:

Components	Amount					Storage
	100 U	500 U	1000 U	5000 U	10000 U	
DNase I	100 U	500 U	1000 U	5000 U	10000 U	-20 °C
DNase I 10× Reaction Buffer	1 ml*1 vial	1 ml*1 vial	1 ml*2 vials	6 ml*1 vial	6 ml*2 vials	-20 °C

Typical protocol for removal of template DNA after *in vitro* transcription


- Add 2 U of DNase I (E00053) to the transcription system for 0.2-1 µg DNA.
- Mix thoroughly, incubate at 37 °C for 30 min.
- Inactivate DNase I by phenol/chloroform extraction.

Typical protocol for removal of genomic DNA from RNA sample


Assemble the reaction in a nuclease-free microcentrifuge tube or PCR strip tube on ice with the following order:

Components	Volume
DNase I (E00053)	10 U
DNase I 10 × Reaction Buffer	10 µl (1 ×)
RNA	~10 µg
RNase free water	to 100 µl
Mix thoroughly, incubate at 37 °C for 30 min.	
Add EDTA with a final concentration of 5 mM and heat the solution at 65 °C for 10 min to inactivate DNase I. Alternatively, DNase I can be inactivated by phenol/chloroform extraction.	
Use the prepared RNA as the template for RT-PCR.	

Data Images



Incubate with 100 ng plasmid DNA at 37 °C for 30 min, tiny amount of GenScript DNase I can completely digest the DNA.



Lane 1: IVT RNA (untreated)
Lane 2: IVT RNA + 2 U DNase I

Add 2 U of DNase I ([E00053](#)) to 20 μ l *in vitro* transcription system, incubate at 37 °C for 30 min, GenScript DNase I can completely digest the DNA template.

References

1. Kienzle, N., et al. "DNase I treatment is a prerequisite for the amplification of cDNA from episomal-based genes." *BioTechniques* 20.4 (1996): 612-616.
2. Anderson, Stephen. "Shotgun DNA sequencing using cloned DNase I-generated fragments." *Nucleic acids research* 9.13 (1981): 3015-3027.
3. Green, Michael R., Tom Maniatis, and D. A. Melton. "Human β -globin pre-mRNA synthesized *in vitro* is accurately spliced in *Xenopus* oocyte nuclei." *Cell* 32.3 (1983): 681-694.

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.