

Version: 03

## DATASHEET

Update: 01/05/2024

# GenCRISPR™ Cas13a (C2c2) Nuclease

Cat. No.: Z03486

### Product Introduction

GenCRISPR™ Cas13a (C2c2) Nuclease is an RNA-guided RNA endonuclease from *Leptotrichia wadei*. This recombinant nuclease is purified from *E. coli* expression and contains C-terminal His Tag. Cas13a (previously known as C2c2) belongs to the Class 2 Type VI CRISPR-Cas system and contains a single protein effector with two HEPN domains. Cas13a exhibits two distinct ribonuclease activities. Cas13a possesses sequence-specific RNA cleavage activity *in vitro* and *in vivo*. And after specific cleavage, Cas13a activates collateral cleavage activity towards adjacent non-specific ssRNA sequences. The two distinct ribonuclease activities of Cas13a enable guide-RNA processing and RNA detection. The RNA detection empowered with Cas13a-crRNA complex and quenched fluorescent ssRNA as reporter has been reported to be versatile, rapid and high-specific.

**Source:** Recombinant Cas13a expressed by *E. coli*

**Species:** *Leptotrichia wadei*

**Accession#:** U2PSH1

**Tag:** C-terminal His Tag

**Apparent Molecular Weight:** ~140 kDa, on SDS-PAGE under non-reducing conditions

**Concentration:** Please refer to the COA for the

specific lot.

**Storage Buffer:** 10 mM Tris-HCl, 300 mM NaCl, 1.0 mM DTT, 50% glycerol, pH 8.0.

**Storage & Stability:** Store at -20 °C for up to 12 months from the date of manufacture. Avoid repeated freeze-thaw cycles. **Do not store below -20 °C!**

**Application:** crRNA-dependent RNA cleavage

### Quality Control Specifications

Assay	Specifications
Appearance	Clear, colorless liquid
Purity by SDS-PAGE	≥ 95%
Activity	≥ 85% by <i>in vitro</i> digestion of ssRNA
Concentration	4 mg/ml±10% as analyzed by A280
Residual DNase	Undetectable
Residual RNase	Undetectable
Endotoxin Level	< 0.2 EU/μg

## Reagents Supplied

Components	Size			Storage
	100 µg	500 µg	1 mg	
GenCRISPR™ Cas13a (C2c2) Nuclease	100 µg	500 µg	1 mg	-20 °C
10 × Cas13a Reaction Buffer*	1.5 ml	1.5 ml	1.5 ml	-20 °C

\* The reaction buffer is intended for Cas13a/crRNA-based Fluorescent Assays described in this manual, and not intended for other applications.

## Protocols for Cas13a/crRNA-Based Fluorescent Assays

### 1. Reagents preparations

- 1) Prepare the specific crRNA, we recommend CRISPR Synthetic Guide RNA Services ([see www.genscript.com](http://www.genscript.com)) for crRNA synthesis.
- 2) Synthesize specific single-strand RNA (ssRNA) as substrate for cleavage; we recommend custom RNA synthesis using GenScript's RNA Synthesis Service ([see www.genscript.com](http://www.genscript.com)).
- 3) Synthesize quenched fluorescent ssRNA; we recommend custom RNA synthesis using GenScript RNA Synthesis Service ([see www.genscript.com](http://www.genscript.com)).
- 4) Prepare the crRNA work solution at an appropriate concentration (e.g. 7.5 ng/µl) by diluting the stock solution with DEPC-Treated Water (nuclease-free) on ice.
- 5) Prepare the substrate ssRNA work solution at an appropriate concentration (e.g. 5 ng/µl) by diluting the stock solution with DEPC-Treated Water (nuclease-free) water on ice.
- 6) The reaction volume is recommended to set at 50 µl, and the reaction conditions should be optimized according to customers' specific applications.

### 2. Assay procedures

- 1) Assemble the reaction in a nuclease-free PCR strip tube for step 1 at room temperature on a clean bench as the following order in the table.
- 2) Assemble the reaction in a black ELISA Plate for step 2 at room temperature on a clean bench in the order indicated in the table.

Step 1: Prepare the Cas13a/crRNA complex	
Components	Volume
10 × Cas13a Reaction Buffer	1 µl
15 ng crRNA	2 µl (7.5 ng/µl)
200 ng GenCRISPR™ Cas13a (C2c2) Nuclease (Z03486)	2 µl (100 ng/µl)
DEPC-Treated Water (nuclease-free)	5 µl
Mix thoroughly, incubate for 10 min at 37 °C in a PCR instrument.	
Step 2: Fluorescent Assay	
Components	Volume
10 ng ssRNA	2 µl (5 ng/µl)
10 × Cas13a Reaction Buffer	4 µl
DEPC-Treated Water (nuclease-free)	29 µl
Cas13a/crRNA complex (from step 1)	10 µl

10 pmol Fluorescent RNA reporter	5 $\mu$ l
Mix thoroughly, incubate for 30 min at 37 °C in an incubator. After incubation, read the fluorescence intensity (Excitation wavelength: 494 nm; Emission wavelength: 518 nm) by a Microplate Reader.	

**References:**

1. Abudayyeh, Omar O., et al. "C2c2 is a single-component programmable RNA-guided RNA-targeting CRISPR effector." *Science* 353.6299 (2016).
2. Gootenberg, Jonathan S., et al. "Nucleic acid detection with CRISPR-Cas13a/C2c2." *Science* 356.6336 (2017): 438-442.
3. Gootenberg, Jonathan S., et al. "Multiplexed and portable nucleic acid detection platform with Cas13, Cas12a, and Csm6." *Science* 360.6387 (2018): 439-444.
4. Iwasaki, Roman S., and Robert T. Batey. "SPRINT: a Cas13a-based platform for detection of small molecules." *Nucleic acids research* 48.17 (2020): e101-e101.

**For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.**

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28 Yongxi Road, Jiangning District, Nanjing, Jiangsu, China