

Version: 03
 Update: 08/24/2023

DATASHEET

Benz-Neburase™, tag-free


Cat. No.: Z03695-10; Z03695-100; Z03695-500

Size: 10 kU; 100 kU; 500 kU


Product Introduction

Species:	<i>Serratia Marcescens</i>
Tag:	Tag-free
Purity:	≥ 95% as analyzed by reducing SDS-PAGE
Endotoxin Level:	≤ 0.1 EU/kU as determined by gel clotting method
Specific Activity:	≥ 1.1×10 ⁶ U/mg One unit of Benz-Neburase, tag-free is defined as the amount of enzyme for a ΔA260 of 1.0 (equivalent to the complete digestion of 37 μg DNA) in 30 min.
Enzyme Activity:	≥ 250 U/μL
Expression System:	<i>E. coli</i>
Theoretical Molecular Weight:	27.5 kDa
Apparent Molecular Weight:	~27.5 kDa, on SDS-PAGE under reducing conditions
Application:	The Benz-Neburase™ is an endonuclease capable of removing all forms of DNA and RNA, including double stranded, single stranded, linearized, and circular forms. The Benz-Neburase is commonly used in biopharmaceutical production such as vaccine, viral vector, gene and cell therapy manufacturing facilities.
Application Note:	The activity of Benz-Neburase™, tag-free requires 1-2 mM Mg ²⁺ .
Formulation:	Supplied as a solution of 20 mM Tris-HCl, 2 mM MgCl ₂ , 20 mM NaCl, 50% Glycerol, pH 8.0.
Storage & Stability:	This product remains stable for up to 2 weeks at 4 °C or up to 24 months at -20 °C. Avoid repeated freeze-thaw cycles. Do not store below -20 °C!

Data Images



Lane M: DNA marker
 Lane 1: PCR product
 Lane 2: GenScript Benz-Neburase™, tag-free + PCR product
 Lane 3: Competitor endonuclease+ PCR product
 Lane 4: Genomic DNA
 Lane 5: GenScript Benz-Neburase™, tag-free + Genomic DNA
 Lane 6: Competitor endonuclease + Genomic DNA
 Lane 7: Plasmid DNA
 Lane 8: GenScript Benz-Neburase™, tag-free + Plasmid DNA
 Lane 9: Competitor endonuclease + Plasmid DNA
 Lane 10: RNA
 Lane 11: GenScript Benz-Neburase™, tag-free + RNA
 Lane 12: Competitor endonuclease + RNA



Lane M: SDS-PAGE marker

Lane R: Reducing (R)

Lane N: Non-reducing (NR)

Purity: > 95% as analyzed by SDS-PAGE

Background

Target Background: The Benz-Neburase is a genetically engineered endonuclease from *Serratia marcescens* used as a DNA eraser in the purification processes of biological molecules. The enzyme cleaves all forms of DNA and RNA into smaller nucleotides of around 5-8 base pairs. Benz-Neburase requires divalent cation, preferably Mg^{2+} for activity, displays a broad pH tolerance, ranging from pH 6 to pH 10, with an optimal pH of 8-8.5, and has a wide temperature tolerance, ranging from 35 °C to 44 °C. The nuclease is a physiologic homodimer and functions more progressively than the monomer. Two disulfide bonds in the nuclease are crucial to its activity and stability. The enzyme is active in a broad range of conditions and is free of proteolytic activity. This makes the enzyme especially useful for biopharmaceutical applications with contaminating DNA residue, such as lysed host cells in viral vector manufacturing processes.

Synonyms: alternative to Benzonase®

Benzonase is a registered trademark of Merck KGaA

References:

1. Nestle, Marion, and W. K. Roberts. "An extracellular nuclease from *Serratia marcescens*: I. Purification and some properties of the enzyme." *Journal of Biological Chemistry* 244.19 (1969): 5213-5218.
2. Benedik, Michael J., and Ulrich Strych. "Serratia marcescens and its extracellular nuclease." *FEMS microbiology letters* 165.1 (1998): 1-13.
3. Filimonova, Maria N., Kurt L. Krause, and Michael J. Benedik. "Kinetic studies of the *Serratia marcescens* extracellular nuclease isoforms." *Biochemistry and molecular biology international* 33.6 (1994): 122-1032.
4. Friedhoff, Peter, et al. "A procedure for renaturation and purification of the extracellular *Serratia marcescens* nuclease from genetically engineered *Escherichia coli*." *Protein expression and purification* 5.1 (1994): 37-43.
5. Franke, Ingo, Gregor Meiss, and Alfred Pingoud. "On the Advantage of Being a Dimer, a Case Study Using the Dimeric *Serratia* Nuclease and the Monomeric Nuclease from *Anabaena* sp. Strain PCC 7120." *Journal of Biological Chemistry* 274.2 (1999): 825-832.

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

生产商: 南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路 28 号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28 Yongxi Road, Jiangning District, Nanjing, Jiangsu, China

GenScript USA, Inc.

860 Centennial Ave. Piscataway, NJ 08854

Tel: 1-732-885-9188