

Version: 02

Update: 02/18/2024

DATASHEET**GenCRISPR™ LbCas12a Nuclease**

Cat. No.: Z03753

Product Introduction

GenCRISPR™ LbCas12a Nuclease is an RNA-guided DNA endonuclease from *Lachnospiraceae bacterium*. Cas12a (previously known as Cpf1) belongs to the Class 2 Type V CRISPR/Cas system. Different from CRISPR/Cas9 system, the ribonucleoprotein (RNP) complex of CRISPR/Cas12a system is formed by Cas12a and crRNA, it recognizes a T-rich protospacer adjacent motif (PAM) and results in a staggered DNA double-strand break (DSB). After the specific cleavage, Cas12a can also activate collateral cleavage activity towards adjacent non-specific ssDNA sequences. Hence, Cas12a nuclease is a good alternative for Cas9 in certain target DNA editing, and provides a novel strategy for DNA detection.

Source: Recombinant Cas12a with a C-terminal NLS expressed by *E.coli*

Species: *Lachnospiraceae bacterium*

Accession#: 6KL9_A

Tag: C-terminal 6× His Tag

Apparent Molecular Weight: 146 kDa, on SDS-PAGE under non-reducing conditions

Concentration: 4 mg/ml

Storage Buffer: 10 mM Tris-HCl, 300 mM NaCl, 0.5 mM DTT, 50% glycerol, pH 7.4.

Storage & Stability: Store at -20 °C for up to 12 months from the date of manufacture. Avoid repeated freeze-thaw cycles. **Do not store below -20 °C!**

Application: crRNA-dependent DNA cleavage

Quality Control Specifications

Assay	Specifications
Appearance	Clear, colorless liquid
Purity by SDS-PAGE	≥ 90%
Concentration by A280	4(± 10%) mg/ml
Bioactivity (<i>in vitro</i>)	≥ 90%
Residual DNase	Undetectable
Residual RNase	Undetectable
Endotoxin Level	< 0.1 EU/μg

Reagents Supplied

Components	Size			Storage
GenCRISPR™ LbCas12a Nuclease	100 µg	500 µg	1 mg	-20 °C
10 × Cas12a Reaction Buffer*	1.5 ml	1.5 ml	1.5 ml	-20 °C

* The reaction buffer is intended for *In vitro* digestion assays described in this manual only, and not intended for other applications.

Protocols for *In vitro* digestion of DNA

1. Reagents preparations

- 1) Synthesize specific crRNA; we recommend GenScript's CRISPR Synthetic Guide RNA Services (see www.genscript.com) for crRNA synthesis.
- 2) Synthesize specific substrate DNA for cleavage; we recommend custom DNA synthesis using GenScript's DNA Synthesis Service (see www.genscript.com), or customers can also apply PCR method to obtain substrate DNA (**Note:** the PCR product needs to be purified prior to further assays).
- 3) Refer to the following procedure to anneal crRNA before preparing the crRNA working solution.

95 °C	5 min	-
94 °C	10s	Decrease 1 °C per cycle, and 69 cycles in total
25 °C		
4 °C	Hold	-

- 4) Prepare the crRNA working solution at an appropriate concentration (e.g., 5 ng/µl) by diluting the stock solution with DEPC-Treated Water (nuclease-free) on ice.
- 5) Prepare the substrate DNA working solution at an appropriate concentration (e.g., 30 ng/µl) by diluting the stock solution with DEPC-Treated Water (nuclease-free) on ice.
- 6) The reaction volume is recommended to be set at 20 µl, and the reaction conditions should be optimized according to your specific applications.

2. Assay procedures

- 1) Assemble the reaction in a nuclease-free PCR strip tube at room temperature on a clean bench, mix the following reagents thoroughly, and incubate for 10 min at 25 °C.

Components	Volume
10 × Cas12a Reaction Buffer	2 µl
10 ng crRNA	2 µl (5 ng/µl)
100 ng GenCRISPR™ LbCas12a Nuclease (Cat. No. Z03753)	2 µl (50 ng/µl)
DEPC-Treated Water (nuclease-free)	12 µl

- 2) Add substrate DNA into the mixture and mix again thoroughly, then incubate for 30 min at 37 °C;

Components	Volume
60 ng Substrate DNA	2 µl (30 ng/µl)
Total reaction volume	20 µl

- 3) Add 1 µl Proteinase K (20 µg/µl), incubate for 30 min at 55 °C. After incubation, determine the digestion efficiency by agarose gel electrophoresis.

Product Validation




Figure 1: *In vitro* digestion efficiency analysis by agarose gel electrophoresis. A 20 μ l reaction in 1 \times Cas12a Nuclease Reaction Buffer containing 60 ng linearized plasmid, 10 ng crRNA, and 100 ng GenCRISPR™ LbCas12a Nuclease (Cat. No. Z03753) for 30 min at 37°C results in a digestion efficiency of linearized plasmid higher than 90%, as determined by agarose gel electrophoresis.

References

1. Zetsche, Bernd, *et al.* "Cpf1 is a single RNA-guided endonuclease of a class 2 CRISPR-Cas system." *Cell* 163.3 (2015): 759-771.
2. Ledford, Heidi. "Alternative CRISPR system could improve genome editing." *Nature News* 526.7571 (2015): 17.
3. Chen, Janice S., *et al.* "CRISPR-Cas12a target binding unleashes indiscriminate single-stranded DNase activity." *Science* 360.6387 (2018): 436-439.

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

生产商：南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路 28 号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28 Yongxi Road, Jiangning District, Nanjing, Jiangsu, China

GenScript USA, Inc.

860 Centennial Ave. Piscataway, NJ 08854

Tel: 1-732-885-9188