


Positive

MS Spectrum



Sample Information		Interface		Prerod Bias	
User	:xu	Nebulizing Gas Flow	:ESI	Detector	:+4.5kv
Date Acquired	: 2008-3-10 17:46:38	CDL Temp	:1.50 L/min	T.Flow	:1.5kv
Sample Name	:PA(224-233),Influenza	CDL Volt	: 250°C	B.conc	: 0.2ml/min
Cat No.	:RP19991	Block Temp	: -25.0V		: 50%H2O/50%MeOH
Injection Volum	: 1.000μL		: 200°C		
Theoretical MW	:1185.30				
Observed MW	:1184.45				