


Positive

MS Spectrum



Sample Information		Interface		Prerod Bias	
User	:Grace	Nebulizing Gas Flow	:1.50 L/min	Detector	:1.5kv
Date Acquired	: 02/15/2009 14:01:58	CDL Temp	: 250°C	T.Flow	: 0.2ml/min
Sample Name	:Beta-Amyloid(1-9)	CDL Volt	: -25.0V	B.conc	: 50%H2O/50%MeOH
Cat No.	: RP20160	Block Temp	: 200°C		
Lot No.	: CA05902				
Injection Volum	: 1.000µL				
Theoretical MW	: 1033.00				
Observed MW	: 1032.60				