


Positive

MS Spectrum



Sample Information		Interface		Prerod Bias	:+4.5kv
User	:Grace	Nebulizing Gas Flow	:1.50 L/min	Detector	:1.5kv
Date Acquired	: 02/23/2009 11:41:25	CDL Temp	: 250°C	T.Flow	: 0.2ml/min
Sample Name	:Beta-Amyloid Protein Precursor 770(135-155)	CDL Volt	: -25.0V	B.conc	: 50%H2O/50%MeOH
Cat No.	: RP20196	Block Temp	: 200°C		
Lot No.	: CA08902				
Injection Volum	: 1.000µL				
Theoretical MW	: 2617.01				
Observed MW	: 2618.25				