

## CERTIFICATE OF ANALYSIS

### Product information

---

Product Name	THE™ Anti-c-Myc [HRP], mAb, Mouse
Product Cat. No.	A00863
Lot No.	10H000181
Quantity	Specified on product label
Presentation	Lyophilized powder from antibody solution with PBS, pH 7.4, containing 1% BSA and 0.01% thimerosal.
Conjugation	Horseradish peroxidase (HRP)
Storage	The product should be stored lyophilized until use. It remains stable in lyophilized form for 3 years from the tested date 08/27/2010. The reconstituted antibody can be stored for 2-3 weeks at 2-8°C or for up to 12 months at -20°C or below. Avoid repeated freezing and thawing cycles.

---

### Quality Control

---

Test Items	Specifications	Results
Reactivity	ELISA titer $\geq$ 1:20,000 Specific to c-Myc tagged fusion protein by WB	Passed

---

### Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:   
Antibody Production Director      Date: 08/27/2010