


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 08/29/16
Time Processed : 14:08:18
Injection Volume : 0.2
Sample Name : Competence-Stimulating Peptide-12261
Cat No. : RP20427
Lot No. : PE4661608
Theoretical MW : 2153.50
Observed MW : 2153.10

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250 °C
Block Temp : 200 °C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH