


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 05/14/16
Time Processed : 14:00:07
Injection Volume : 0.2
Sample Name : Biotin-TAT (47-57)
Cat No. : RP20408
Lot No. : PE4351605
Theoretical MW : 1786.20
Observed MW : 1786.50

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH