


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 04/20/16
Time Processed : 14:18:08
Injection Volume : 0.2
Sample Name : Dermorphin Analog
Cat No. : RPI3522
Lot No. : PE7321604
Theoretical MW : 901.43
Observed MW : 901.60

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH