


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 02/26/16  
Time Processed : 15:24:47  
Injection Volume : 0.5  
Sample Name : Myelin Basic Protein (MBP) (synthetic)  
Cat No. : RPI0314  
Lot No. : PE1941602  
Theoretical MW : 1390.50  
Observed MW : 1390.80

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH