


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 08/28/16  
Time Processed : 17:16:07  
Injection Volume : 0.2  
Sample Name : Competence-Stimulating Peptide-12261  
Cat No. : RP20427  
Lot No. : PE9891608  
Theoretical MW : 2153.50  
Observed MW : 2152.80

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH