


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/04/16
Time Processed : 13:28:30
Injection Volume : 0.1
Sample Name : β -Amyloid (1-42), human
Cat No. : RP10017
Lot No. : PE9001605
Theoretical MW : 4514.10
Observed MW : 4513.70

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250 °C
Block Temp : 200 °C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH