

Version: 02

Update: 12/13/2022

CERTIFICATE OF ANALYSIS

Product Information

Product Name:	CytoSinct™ CD3 Nanobeads, human
Cat. No.:	L00896-0.5
Lot No.:	R100324ABB
Storage:	Store at 2-8 °C
Expiry Date:	10/18/2024

Quality Control

Test Items	Specifications	Results
Physical appearance	Amber to brown liquid with no visible impurities	Pass
Selection purity	Satisfactory selection purity from PBMC	Pass
Selection recovery	Satisfactory selection recovery from PBMC	Pass

Date: 04/17/2024

Approved by: 存燕

Quality Department

For research use only. Not intended for human and animal therapeutic or diagnostic use.

生产商: 南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路 28 号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28 Yongxi Road, Jiangning District, Nanjing, Jiangsu, China