

Version: 02

Update: 12/13/2022

CERTIFICATE OF ANALYSIS

Product Information

| | |
|---------------|---------------------------------|
| Product Name: | CytoSinct™ CD3 Nanobeads, human |
| Cat. No.: | L00896-0.5 |
| Lot No.: | R100324BBB |
| Storage: | Store at 2-8 °C |
| Expiry Date: | 11/26/2024 |

Quality Control

| Test Items | Specifications | Results |
|---------------------|--|---------|
| Physical appearance | Amber to brown liquid with no visible impurities | Pass |
| Selection purity | Satisfactory selection purity from PBMC | Pass |
| Selection recovery | Satisfactory selection recovery from PBMC | Pass |

Date: 04/17/2024

Approved by: 陈燕

Quality Department

For research use only. Not intended for human and animal therapeutic or diagnostic use.

生产商: 南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路 28 号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28 Yongxi Road, Jiangning District, Nanjing, Jiangsu, China