


Mass Spectrum



Sample Information

Month-Day Processed : 11/15/20
Time Processed : 01:18:59 PM
Injection Volume : 0.4
Sample Name : E1, (E257)
Cat No. : RP20441
Lot No. : E4632820K
Theoretical MW : 1005.17
Observed MW : 1004.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH