


# Mass Spectrum



## Sample Information

User : Gary  
 Month-Day Processed : 08/23/18  
 Time Processed : 2:02:22 PM  
 Injection Volume : 0.8  
 Sample Name : ANP (1-28), human, porcine  
 Cat No. : RP11927  
 Lot No. : PE5401808  
 Theoretical MW : 3080.50  
 Observed MW : 3080.1

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250|æ  
 Block Temp : 200|æ

Interface Bias : +4.5kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2ml/min  
 B.conc : 50%H2O/50%MeOH