


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 01/13/17
Time Processed : 2:24:39
Injection Volume : 1
Sample Name : HA Peptide
Cat No. : RP11735
Lot No. : PE7201701
Theoretical MW : 1102.15
Observed MW : 1101.55

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH