


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	02/20/19	Interface Bias	: +4.5kV
Time Processed :	08:38:48 PM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.2	T.Flow	: 0.2ml/min
Sample Name :	Bradykinin	B.conc	: 50%H2O/50%MeOH
Cat No. :	RP11129		
Lot No. :	PE3451902		
Theoretical MW :	1060.20		
Observed MW :	1060.6		