


Mass Spectrum



Sample Information
Month-Day Processed : 08/26/21
Time Processed : 23:24:54
Injection Volume : 0.2
Sample Name : Secretin, human
Cat. No. : RP10249
Lot. No. : E9466721H
Theoretical MW : 3039.50
Observed MW : 3039.6

Interface : ESI
Nebulizing Gas Flow : 1.5 L/min
CDL Temp : 250
Block Temp : 200
Equipment : ZJ21010035
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH