


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 12/14/16
Time Processed : 11:00:08
Injection Volume : 1
Sample Name : Melitin
Cat No. : RPI0290
Lot No. : PE5501612
Theoretical MW : 2846.46
Observed MW : 2846.60

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH