


Mass Spectrum



Sample Information

Acquired by : Gary  
Month-Day Processed : 12/02/17  
Time Processed : 21:27:47  
Injection Volume : 0.2  
Sample Name : Eptifibatide  
Sample ID : U0297CJ080-1  
Theoretical MW : 832.01  
Observed MW : 831.50

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH