


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 06/02/17  
Time Processed : 11:38:29  
Injection Volume : 0.2  
Sample Name :  $\beta$ -Amyloid (1-40)  
Cat No. : RP10004  
Lot No. : PE2741705  
Theoretical MW : 4329.82  
Observed MW : 4330.05

## Interface

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

## Interface Bias

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH