


Mass Spectrum



Sample Information

Month-Day Processed : 08/10/21
Time Processed : 20:46:00
Injection Volume : 0.4
Sample Name : Somatostatin
Cat No. : RP10226
Lot No. : E4117321H
Theoretical MW : 1637.88
Observed MW : 1637.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH