


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 06/04/18
Time Processed : 15:49:38
Injection Volume : 0.2
Sample Name : HIV - 1 Rev (34 - 50)
Cat No. : RP20362
Lot No. : PE6071805
Theoretical MW : 2437.75
Observed MW : 2437.95

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH