


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 07/30/18
Time Processed : 18:56:02
Injection Volume : 0.2
Sample Name : Hexa-His
Cat. No. : RP11737
Lot. No. : PE1761807
Theoretical MW : 840.85
Observed MW : 840.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH