


Mass Spectrum



Sample Information

Acquired by: Gary
Month-Day Processed: 07/17/19
Time Processed: 02:41:18 PM
Injection Volume: 0.3
Sample Name: Arg-Gly-Asp-Cys
Cat. No.: RP10860
Lot. No.: E5457819G
Theoretical MW: 449.48
Observed MW: 449.2

Interface: ESI
Nebulizing Gas Flow: 1.5L/min
CDL Temp: 250
Block Temp: 200

Equipment: GK11010007
Interface Bias: +4.5 kV
Drying Gas Flow: 5 L/min
T.Flow: 0.2 ml/min
B.conc: 50%H2O/50%MeOH