


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 09/20/19  
Time Processed : 09:45:38 AM  
Injection Volume : 0.4  
Sample Name : Protease-Activated Receptor-4  
Cat No. : RP20296  
Lot No. : E4865519I  
Theoretical MW : 666.77  
Observed MW : 666.8

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200