


# Mass Spectrum



## Sample Information

User : Gary  
 Month-Day Processed : 12/29/16  
 Time Processed : 0:02:38  
 Injection Volume : 1  
 Sample Name : Brain Natriuretic Peptide (BNP) (1-32), rat  
 Cat No. : RPI1121  
 Lot No. : PE3441612  
 Theoretical MW : 3452.94  
 Observed MW : 3453.00

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250°C  
 Block Temp : 200°C

Interface Bias : +4.5kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2ml/min  
 B.conc : 50%H2O/50%MeOH