


Mass Spectrum



Sample Information
Month-Day Processed : 05/09/21
Time Processed : 10:10:55 AM
Injection Volume : 0.4
Sample Name : Hexa-His
Cat. No. : RP11737
Lot. No. : E2666321E
Theoretical MW : 840.85
Observed MW : 840.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH