


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 05/28/18
 Time Processed : 8:30:33
 Injection Volume : 0.2
 Sample Name : β -Amyloid (1-40)
 Cat. No. : RP10004
 Lot. No.: PE1381804
 Theoretical MW : 4329.82
 Observed MW : 4329.75

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H_2O /50%MeOH