


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 01/21/19  
Time Processed : 08:23:13 PM  
Injection Volume : 0.2  
Sample Name : Arg-Gly-Asp-Ser  
Cat No. : RP10861  
Lot No. : PE8861901  
Theoretical MW : 433.42  
Observed MW : 433.2

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH