


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 07/02/19
Time Processed : 07:09:58 PM
Injection Volume : 0.3
Sample Name : Calcitonin Gene Related Peptide (CGRP) (8-37), human
Cat No. : RP11089
Lot No. : E0909519G
Theoretical MW : 3125.59
Observed MW : 3126.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH