


Mass Spectrum



Sample Information

Month-Day Processed : 08/12/20
Time Processed : 08:09:23 AM
Injection Volume : 0.3
Sample Name : Dynorphin A(1-13), porcine
Cat No. : RP20492
Lot No. : E0583720H
Theoretical MW : 1604.00
Observed MW : 1603.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH