


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 08/01/19
Time Processed : 06:50:10 PM
Injection Volume : 0.3
Sample Name : Melanotan (MT)-II
Cat No. : RP17323
Lot No. : E5643519H
Theoretical MW : 1024.18
Observed MW : 1024.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH