


Mass Spectrum



Sample Information
Acquired by : Gary
Month-Day Processed : 11/01/17
Time Processed : 22:47:32
Injection Volume : 0.2
Sample Name : Urocortin II, human
Sample ID : U1292CI200-1
Theoretical MW : 4450.28
Observed MW : 4451.10

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH

: +4.5 kV