


# Mass Spectrum



Sample Information  
User : Gary  
Month-Day Processed : 12/01/19  
Time Processed : 12:43:07 PM  
Injection Volume : 0.4  
Sample Name : Galanin, human  
Cat No. : RP10801  
Lot No. : E7321319K  
Theoretical MW : 3157.41  
Observed MW : 3156.5

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH