


# Mass Spectrum



Sample Information

User : Gary

Month-Day Processed : 09/19/19

Time Processed : 09:23:57 AM

Injection Volume : 0.3

Sample Name : Secretoneurin, rat

Cat No. : RP13702

Lot No. : E4908019I

Theoretical MW : 3651.95

Observed MW : 3651.8

Equipment : GK11010007

Interface Bias : +4.5kV

Drying Gas Flow : 5 L/min

T.Flow : 0.2ml/min

B.conc : 50%H2O/50%MeOH

Interface : ESI

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250

Block Temp : 200