


# Mass Spectrum



Sample Information  
 Acquired by: Gary  
 Month-Day Processed: 04/10/17  
 Time Processed: 5:31:59  
 Injection Volume: 1  
 Sample Name: Somatostatin  
 Cat No.: RPI0226  
 Lot No.: PE4961703  
 Theoretical MW: 1637.88  
 Observed MW: 1638.30

Interface :ESI  
 Nebulizing Gas Flow :1.5L/min  
 CDL Temp :250°C  
 Block Temp :200°C

Interface Bias :+4.5 kV  
 Drying Gas Flow :5 L/min  
 T.Flow :0.2 ml/min  
 B.conc :50%H2O/50%MeOH