


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/18/19
Time Processed : 02:59:59 PM
Injection Volume : 0.3
Sample Name : PACAP (6-38), human, ovine, rat
Cat No. : RP10338
Lot No. : E4700119I
Theoretical MW : 4024.76
Observed MW : 4024.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH