


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 10/21/18  
Time Processed : 17:53:31  
Injection Volume : 0.2  
Sample Name : Myelin Basic Protein (MBP) (synthetic)  
Cat. No. : RPI0314  
Lot. No. : PE8021810  
Theoretical MW : 1390.50  
Observed MW : 1390.8

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH