


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 01/03/17  
Time Processed : 11:25:38  
Injection Volume : 1  
Sample Name :  $\beta$ -Amyloid (1-42), human  
Cat No. : RP10017  
Lot No. : PE9781611  
Theoretical MW : 4514.10  
Observed MW : 4514.00

## Interface

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

## Interface Bias

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H<sub>2</sub>O/50%MeOH