


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/21/19
Time Processed : 09:53:49 AM
Injection Volume : 0.3
Sample Name : Neuromedin N
Cat No. : RP10509
Lot No. : E4747119I
Theoretical MW : 745.95
Observed MW : 745.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH