


# Mass Spectrum



Sample Information  
Month-Day Processed : 11/01/20  
Time Processed : 05:50:09 PM  
Injection Volume : 0.4  
Sample Name : C-Peptide (3-33), human  
Cat. No. : RP11116  
Lot. No. : E9806820J  
Theoretical MW : 3020.26  
Observed MW : 3020.1

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH