


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/03/19
Time Processed : 11:37:05 AM
Injection Volume : 0.2
Sample Name : Peptide YY (PYY) (3-36), human
Cat No. : RP10354
Lot No. : PE8151902
Theoretical MW : 4049.48
Observed MW : 4049.0

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH