


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 05/26/19
Time Processed : 07:42:08 PM
Injection Volume : 0.3
Sample Name : Hsp70-derived octapeptide
Cat No. : RP20220
Lot No. : PE8551905
Theoretical MW : 858.91
Observed MW : 858.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH