


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 11/01/18  
Time Processed : 16:43:41  
Injection Volume : 0.2  
Sample Name : Neuropeptide Y, human  
Cat No. : RP10492  
Lot No. : PE5671810  
Theoretical MW : 4271.68  
Observed MW : 4271.4

## Interface

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

## Equipment

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH